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MASAŻERA

KARKU

MODEL: ORO-NECK MASSAGER

SPIS TREŚCI

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USAGE CONDITION

1. Store in dry,dust free location .
2. Do not spill liquid on the device.
3. Do not immerse the device in water.
4. Do not clean with chemicals.

SAFETY PRECAUTIONS

Presons with the following condition must consult the doctor before using this unit:

1. Acute disease.
2. Malignant tumor.
3. Infectious disease.
4. Pregnancy.
5. Cardiac dusfunction.
6. High fever.
7. Abnormal blood pressure.
8. Skin sensory disorders or skin problems.
9. Receiving medical treatment.especially those feeling discomfort.

Should not be used by persons in the first trimester of pregnancy, fitted with a pacemakers or AI CD.or being treated for an existing deep vein thrombosis.

WARNING

If you feel any problems with your body or skin.consult your medical practitioner an follow their instructions.

Do not this unit while sleeping.

Do not use on your broken skin or damaged sensitive nerve damage.

Always starting with the lower level then increaase the level gradually. For safety.driver should not be using while driving.

DESCRIPTION

Innovative 3D kneading massage balls that hit the acupoints perfectly.
Relieve muscle soreness and stiffness.

Promote blood circulation.

Multifunctional: massage every part of body, from neck to toe. Can be used anywhere.in a car.office and home.

THE RIGHT PLACE FOR THE MASSAGE

Foot

Calf

Legs

Waist

Abdomen

Back

Neck

Back

Waist

Abdomen

Hips

Calf

Foot

Arms

Height adjustment
with arms

BUTTON FUNCTION

- Power
- Direction
- Speed
- Heat

POWER

Press to turn on the massager with heat Auto timer will shut it off after 15minutes. Press the Power bouton anytime to turn the unit off.

DIRECTION

When the unit is turned on.the massage roller will automatically change its rolling direction every 60 seconds. Press the button to change and maintain rotaing direction of the massager.

SPEED

3 different speed from low to high, choose the suitable one according to your need.

HEAT

Press to turn on/off the heat function while massage continues.

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